

## Horizon Technology SPE Disk Cross Reference Guide

EPA Method Number/ Sample Type	Atlantic/Pacific SPE Disk or Carbon Cartridge
<b>500 Series/Drinking Water</b>	
506	C18/HC C18
507	C18/HC C18
508.1	C18/HC C18
515	DVB
522	Carbon Cartridge
525.2/525.3	C18/HC C18/DVB
526	DVB
527	DVB
528	DVB
529	DVB
532	C18/HC C18/HLB-M
548	C18/HC C18/HLB-M
550.1	C18/HC C18
553.1	C18/HC C18
554	C18/HC C18
<b>8000 Series/Ground Water</b>	
8015	DVB
8061	C18/HC C18/DVB
8081	C18/HC C18/DVB/HLB-M
8082	C18/HC C18/DVB/HLB-M
8141	DVB
8151	DVB
8270	HLB-M/H With Carbon Cartridge
8270	8270 One Pass
8315	C18/HC C18/DVB/HLB-M
8318	HLB
8330	HLB
<b>600 Series/Waste Water</b>	
608	C18/HC C18/HC C18
625	HLB-M/H
<b>1600 Series/Waste Water</b>	
1613	C18/HC C18/DVB
1614	C18/HC C18/DVB
1657	C18/HC C18/DVB
1664	Pacific/Pacific Premium
1668	C18/HC C18/DVB
1694	HLB-M/L

Sample Type	Atlantic/Pacific SPE Disk or Carbon Cartridge
<b>Known Endocrine Disruptors</b>	
Organochlorine Pesticides	C18 (608) C18/DVB or HLB-M (8081)
PCBs	DVB or HLB -M
Bisphenol A	C18/HC C18
PBDEs	C18/HC C18
Alkyl phenols	C18/HC C18
Dioxins & Furans	C18/HC C18/DVB
Phthalates	C18/HC C18/DVB
Herbicides	C18/HC C18
Hormones / Sex Steroids	HLB -M
PAHs	C-18, DVB or HLB- M
PPCPs	HLB -M or H
<b>Emerging Contaminants</b>	
Pharmaceutical Compounds	HLB-M/L
Endocrine Disruptors	HLBM/L
Carbaryl	C18/HC C18
Microcystins	C18/HC C18
Nonyl Phenol (surfactant)	DVB
Bisphenol A	DVB

<b>Ordering Information</b>	
Atlantic C18	47-2346-02
Atlantic High Capacity C18	47-2346-13
Atlantic DVB	47-2346-06
Atlantic HLB H	47-2346-10
Atlantic HLB M	47-2346-09
Atlantic HLB L	47-2346-08
Atlantic 8270 One Pass	47-2346-11
Carbon Cartridge	49-2620
Pacific - 1664A - 47 mm	1664-47-HT
Pacific- 1664A - 90 mm	1664-90-HT
Pacific Premium 47 mm	1664-47-PHT
Pacific Premium 100 mm	1664-100-PHT



\*Choose the appropriate disk based on total volume and the capacity of compounds per sample.